



The best technical
ability & constant technical
development.
Provide your company
with a valuable
service & consulting

Vision Scientific Laboratory Instruments

For Research and Testing

Incubator

VS-3155i / VS-3250i

- ✓ Having the best condition on test for over room temperature
- ✓ Digital PID Controller
- ✓ Easy to set by Touch Key
- ✓ Auto tuning function for high temperature precision
- ✓ Safety device for over temperature heating
- ✓ Forced Circulation type for temperature distribution with high precision
- ✓ Inner Glass door to be reduced the temperature change
- ✓ 155Liter and 250Liter Capacity
- ✓ 10 Step Temperature Program Control (Option, '-P' Model)

Incubator	VS-3155i	VS-3250i
Dimension (in)	480 × 460 × 700(H)mm	500 × 500 × 1,000(H)mm
(out)	640 × 600 × 955(H)mm	650 × 680 × 1,340(H)mm
Capacity	155 Liter	250 Liter
Weight	65kg	94kg
Temp. Controller	Digital PID Controller	
Temp. Range	Ambient +5°C to 70°C (at 25°C)	
Temp. Accuracy	±0.2°C at 37°C	
Display	Dual LED Display	
Material (in)	Stainless Steel	
(out)	Steel Plate with Powder Coating	
Door	Double Door, Silicone Packing Magnet Door	
Inner Door	Tempered Safety Glass Door	
Circulation Fan	Internal Convection Fan	
Shelves	3EA, Adjustable Type	
Safety Device	Exclusive over Temp. Protector	
Power	220V, 310W, 60Hz	

※ This specification is made basis of 25'C on room temperature

Unit Shape, Specification and Controller can be changed a little without any notice for better quality



VS-3155i

VS-3250i

www.visionsci.co.kr



ISO 9001



ISO 14001



Excellent Machine



VISION SCIENTIFIC CO.,LTD

Head Quarter

(823 Toprip-Dong) 260 Techno 2-Ro, Yuseong-Gu, Daejeon-Si, Korea
TEL.032) 672 – 7779 / FAX. 032) 676-3977

Costechpia 10. Ver. 1.0 on 2015.9.1 VISION SCIENTIFIC